

Certificate

Manufacturer	ResMed Paris 240 rue de la Motte 77550 Moissy-Cramayel France
Device	VS III™ Lung ventilator for medical use

To whom it may concern,

The US Department of Transportation (DOT) Final Rule, "*Nondiscrimination on the Basis of Disability in Air Travel*" (73 FR 27614 which updates Title 14 CFR Part 382), effective May 13, 2009 provides important new requirements for the accommodation of passengers with respiratory assistive devices.

In line with these requirements, respiratory assistive devices may be used onboard an aircraft, without further testing by the carrier, provided they have been tested for Electromagnetic Compatibility (EMC) in accordance with the current version of RTCA/DO-160, Section 21, Category M. ResMed recommends that customers still check with their airline before booking a flight.

ResMed has successfully completed testing for the VS III™ respiratory assistive device.

VS III™ powered by its internal battery or by a ResMed Power Station II (RPS II) complies with RTCA/DO-160G, Section 21, Category M.

Moissy-Cramayel, 22 August, 2011



Jean-Dominique Behety

Operations Director, ResMed Paris